

Development and validation of a modified EUCAST yeast broth microdilution MIC method for rezafungin to mitigate nonspecific binding through incorporation of Tween 20

D. Zuill¹, A. Almaguer¹, J. Donatelli¹, M. C. Arendrup², J. B. Locke^{1*}

¹Cidara Therapeutics, Inc., San Diego, CA, USA; ²Statens Serum Institut, Copenhagen, Denmark

*Jeffrey B. Locke, PhD
jlocke@cidara.com
+1 858 752 6427

INTRODUCTION

- The EUCAST broth microdilution (BMD) methodology generates potent echinocandin MIC values for *Candida* species, which can be problematic at lower drug concentrations if nonspecific binding occurs.
- In a prior multicentre EUCAST study, the novel echinocandin rezafungin (RZF) demonstrated nonspecific binding sufficient to cause significant variability for more susceptible species (e.g. *C. albicans* and *C. tropicalis*) between different MIC plate types.¹
- Herein we investigated the use of surfactants as a means to mitigate nonspecific binding and generate a viable EUCAST BMD method for rezafungin.

METHODS

- Initial antifungal-surfactant range finding checkerboard assays utilized an abbreviated EUCAST BMD "stamping" method (4 μ L surfactant and 4 μ L drug at 50X were dispensed into MIC plates containing 196 μ L of pre-inoculated RPMI-1640 media).
- All other MICs were performed according to EUCAST E.Def 7.3.2.
- Fluconazole was selected as a control for all experiments due to it being minimally impacted by nonspecific binding or surfactants.
- Depending on the experiment, upwards of 3 surfactants (Tween 20 [T20], 80, and Triton X-100), 4 tissue culture-treated MIC plates (Corning-CLS3596, NUNC-167008, Grenier-655180, and Falcon-353072), and a variety of WT and *fks* mutant strains representing 7 *Candida* species (*C. albicans*, *C. glabrata*, *C. tropicalis*, *C. krusei*, *C. parapsilosis*, *C. auris*, and *C. dubliniensis*) were evaluated.

RESULTS

- RZF-surfactant checkerboard assays conducted on a subset of WT and *fks* mutants demonstrated similar abilities of the 3 surfactants to mitigate nonspecific binding (data not shown).
- T20 was selected for further evaluation due to its existing use in EUCAST mould MIC/MEC testing (E.Def 9.3.2) for preparing conidia suspensions.
- 0.002% T20 was the optimal concentration to mitigate RZF nonspecific binding across strains (data not shown).

RESULTS (con't)

Fig. 1 Impact of T20 on WT *Candida* spp. RZF MICs between higher and lower binding plates





Fig. 2 WT and *fks* *C. albicans* and *C. glabrata* RZF MICs in higher and lower binding plates +/- T20



- Comparison of plates exhibiting higher (NUNC) and lower (Falcon) levels of nonspecific binding showed that 0.002% T20 successfully normalized MIC values for panels of WT strains of all 7 *Candida* species (Fig. 1) without compromising the ability to distinguish *fks* mutant and WT strains (Fig. 2).

RESULTS (con't)

Table 1. Pooled RZF MIC distributions across 4 plate types for EUCAST QC strains +/- T20

Strain	T20	MIC (mg/L)															
		0.0005	0.001	0.002	0.004	0.008	0.016	0.031	0.063	0.125	0.3	0.5	1	2	4		
<i>C. albicans</i> ATCC 64548	(-) E1					1	7	13	13	6							
	(+) (+)	2	32	6													
<i>C. albicans</i> ATCC 64550	(-) E1						2	12	19	7							
	(+) (+)		1	27	12												
<i>C. albicans</i> CNM-CL-F8555	(-) E1						1	7	16	16							
	(+) (+)		1	28	11												
<i>C. krusei</i> ATCC 6258	(-) E1								7	23	9	1					
	(+) (+)								7	13							
<i>C. krusei</i> CNM-CL-3803	(-) E1									5	24	11					
	(+) (+)								38	2							
<i>C. parapsilosis</i> ATCC 22059	(-) E1												2	22	16		
	(+) (+)												29	11			

10 independent replicates MICs per strain per plate type; *med*

- T20 normalization of rezafungin MIC values was further validated in a 10-replicate analysis of the 6 EUCAST *Candida* QC strains across all 4 plate types (Table 1).

CONCLUSIONS

- Incorporation of 0.002% T20 into EUCAST BMD MIC assays for rezafungin diminishes nonspecific binding, normalizes MIC values across plate types, and does not impact differentiation of WT vs. *fks* mutant strains.
- The replicate runs of all EUCAST *Candida* spp. QC strains demonstrated high reproducibility of this method across the same plate types used in the prior multicentre study in which bimodal MIC distributions were observed.¹
- This promising methodology adaptation is currently undergoing further validation in a multicentre study and, if successful, could potentially benefit other antifungal agents for which interlaboratory variation, nonspecific binding, and/or plastics choice issues have been documented.^{2,3,4}

REFERENCES

- Arendrup MC, et al. *CMI*, 24 (2016) 1200-1204.
- Arendrup MC, et al. *AAC*, 63 (2019) e00659-19.
- Esping-Ingraff A, et al. *AAC*, 57 (2013) 5836-5842.
- Howard SJ, et al. *AAC*, 57 (2013) 5426-5431.

ACKNOWLEDGMENTS / DISCLOSURES

DZ, AA, JD, JBL: employees and stockholders of Cidara Therapeutics at the time this work was conducted. MCA (outside this study in the past 5 years): research grants/contract work (to affiliation [SS]) from Amplyx, Basilea, Cidara, F2G, Gilead, Novabiotics, Synexis and T20systems; speaker honoraria (personal fee) from Astellas, Chiesi, Gilead, MSD, and SEGES. MCA is the current chairman of the EUCAST-AFST.